

## How to set up an Anaconda environment for DNN pureshift

- Download the GitHub repository
- Open terminal (for Linux/Mac) or Anaconda prompt (for Windows) in the downloaded folder
- Create a conda environment for the environment.yml file; for example, the following command will create an environment named "pureshift"

**conda env create -f environment.yml**

- Then the new environment, **pureshift**, can be activated with the command

**conda activate pureshift**

- In the terminal, type ./run.sh (Linux/MAC) or run.bat in the anaconda prompt (Windows OS), and the output spectra will be stored as a .pdf file, with all individual spectra in the output folder in .csv format.

## How to use DNN pureshift module

**run.sh/ run.bat** looks like the following,

```
python .pureshift \  
    --spin_echo 'NMR/Kanamycin/3/' \  
    --psyche 'NMR/Kanamycin/1002/' \  
    --output 'pred' \  
    --phase 0 \  
    --p0 -125.0 \  
    --p1 355.0 \  
    --pp0 90.4 \  
    --pp1 295.0 \  
    --clear_phase 0 \  
    --scaling 20 \  
    --ver_scale 1.0 \  

```

## Arguments available in the DNN pureshift

Argument	Compulsory / Optional	Purpose	Remarks
--spin_echo	Compulsory	Path of spin-echo spectra	Example: --spin_echo 'NMR/Kanamycin/3/'
--psyche	Optional	Path of PSYCHE spectrum	Example: --psyche 'NMR/Kanamycin/1002/'  <b>if PSYCHE is not available delete this line</b>
--output	Compulsory	The output folder name where the processed spectra will be saved in .csv format	Example: --output 'pred'  It can accept any string variable
--clear_phase	Optional	Resets all phase values to 0.0 in all spectra	Example: --clear_phase 0  Always begin with --clear_phase 1 for each new sample analysis  For any new sample or while performing phase correction, if things go out of control, --clear_phase 1 will reset all phases to 0.0
--phase	Optional	if phase correction is required set --phase = 1, else --phase = 0 (it is also the default and generates pure shift)	Example: --phase 1  Always begin with --phase 1 for every new sample analysis  Once the phase correction is completed, set --phase 0 for DNN pure shift analysis
--p0	Optional	Spin-echo spectra zero-order phase that needs to be changed  It is only used with --phase 1	Example: --p0 125.0  We may not be able to change the zero-order phase by 125.0 in a single step; however, this can be accomplished over multiple steps. For example, consider a phase change of 70.0 in the first step and only 55.0 in the next.
--p1	Optional	Spin-echo spectra first-order phase that needs to be changed  It is only used with --phase 1	Similar to the --p0
--pp0	Optional	PSYCHE zero-order phase that needs to be changed  It is only used with --phase 1	Similar to the --p0
--pp1	Optional	PSYCHE first-order phase that needs to be changed	Similar to the --p0

		It is only used only --phase 1	
--scaling	Optional	Vertical space between the spectra can be changed using this argument	<p>Example: --scaling 20</p> <p>To plot the spectra without any vertical overlap, this argument can be used</p> <p>For example, while plotting Kanamycin spectra, we found that --scaling 20 (also the default setting) clearly separate the spectra</p>
--ver_scale	Optional	This helps in changing the intensity of spectra	<p>Example: --ver_scale 1.0</p> <p>The default value is 1.0 (for DNN pureshift processed spectrum)</p> <p>During plotting, the intensity of DNN pureshift is raised by 1.0, PSYCHE by 64.0, and conventional by 4.0</p>

## How to modify the arguments in run.sh/run.bat file

1. As an initial step, reset all phases to zero with **--clear\_phase 1**



```
python .pureshift \  
    --spin_echo 'NMR/Kanamycin/3/' \  
    --psyche 'NMR/Kanamycin/1002/' \  
    --output 'pred' \  
    --phase 0 \  
    --p0 0.0\  
    --p1 0.0 \  
    --pp0 0.0 \  
    --pp1 0.0 \  
    --clear_phase 1 \  
    --scaling 20 \  
    --ver_scale 1.0 \  

```

2. Along with **--phase 1**, set all of the **p0, p1, pp0, and pp1 to 0.0**, and the spectra will be plotted without any phase correction values. It aids in determining which phase values need to be adjusted. The resulting spectra are shown below.

```
python .pureshift \  
    --spin_echo 'NMR/Kanamycin/3/' \  
    --psyche 'NMR/Kanamycin/1002/' \  
    --output 'pred' \  
    --phase 1 \  
    --p0 0.0\  
    --p1 0.0 \  
    --pp0 0.0 \  
    --pp1 0.0 \  
    --clear_phase 0 \  
    --scaling 20 \  
    --ver_scale 1.0 \  

```





**Spectra after resetting phases to zero with `--clean_phase 1` and plotted with `--phase 1`**

3. Correct the phases for both spin-echo and psyche with **--phase 1**, and set all of the p0, p1, pp0, and pp1 to the expected values. For example, in this scenario, the initial phase values used are **--p0 -75.0, --p1 270.0, --pp0 50.0, and --pp1 180.0**.

```
python .pureshift \  
    --spin_echo 'NMR/Kanamycin/3/' \  
    --psyche 'NMR/Kanamycin/1002/' \  
    --output 'pred' \  
    --phase 1 \  
    --p0 -75.0\  
    --p1 270.0 \  
    --pp0 50.0 \  
    --pp1 180.0 \  
    --clear_phase 0 \  
    --scaling 20 \  
    --ver_scale 1.0 \  

```



**Spectra after the initial step of phase correction with `--phase 1`**



4. Phases to be adjusted further to get proper spectra for both spin-echo and psyche, with **--phase 1**. In this step,
- p0 -50.0 (total phase is  $-50.0 + (\text{from the previous step } -75.0) = -125.0$ )
  - p1 85.0 (total phase is  $85.0 + 270.0 = 355.0$ )
  - pp0 40.4 (total phase is  $40.4 + 50.0 = 90.4$ )
  - pp1 115.0 (total phase is  $115.0 + 180.0 = 295.0$ )

```
python .pureshift \  
    --spin_echo 'NMR/Kanamycin/3/' \  
    --psyche 'NMR/Kanamycin/1002/' \  
    --output 'pred' \  
    --phase 1 \  
    --p0 -50.0\  
    --p1 85.0 \  
    --pp0 40.4 \  
    --pp1 115.0 \  
    --clear_phase 0 \  
    --scaling 20 \  
    --ver_scale 1.0 \  

```

The phase values have now been appropriately set, and this can be seen from the spectra given below. Depending on the requirements, we may need to repeat steps 1–4.






**Spectra after full phase correction**

- Once the phase values have been appropriately adjusted, use the argument **--phase 0**, which will execute the **DNN pure shift processing**. The additional parameters, **--scaling** and **--ver\_scale** (which can be set at default values), will aid in plotting the spectra with appropriate intensities and spacing between them. Finally, the processed spectra in **.csv** format are saved in the output folder along with a **.pdf** file.

```
python .pureshift \
  --spin_echo 'NMR/Kanamycin/3/' \
  --psyche 'NMR/Kanamycin/1002/' \
  --output 'pred' \
  --phase 0 \
  --p0 0.0 \
  --p1 .0 \
  --pp0 0.0 \
  --pp1 0.0 \
  --clear_phase 0 \
  --scaling 20 \
  --ver_scale 1.0 \
```



Final spectra are saved as a PDF file, wherein **--ver\_scale**, and **--scale** arguments are depicted

## How to process only spin-echo spectra while PSYCHE is unavailable

To process only spin-echo spectra with DNN pureshift when PSYCHE is unavailable, delete the `--psyche` argument, which also applies in the phase correction step.

### Phase correction:

```
python .pureshift \  
  --spin_echo 'NMR/Kanamycin/3/' \  
  --psyche 'NMR/Kanamycin/1002/' \ (delete this line)  
  --output 'pred' \  
  --phase 1 \  
  --p0 0.0\  
  --p1 .0 \  
  --pp0 0.0 \  
  --pp1 0.0 \  
  --clear_phase 0 \  
  --scaling 20 \  
  --ver_scale 1.0 \
```

### DNN pureshift processing:

```
python .pureshift \  
  --spin_echo 'NMR/Kanamycin/3/' \  
  --psyche 'NMR/Kanamycin/1002/' \ (delete this line)  
  --output 'pred' \  
  --phase 0 \  
  --p0 0.0\  
  --p1 .0 \  
  --pp0 0.0 \  
  --pp1 0.0 \  
  --clear_phase 0 \  
  --scaling 20 \  
  --ver_scale 1.0 \
```

## How to acquire spin-echo NMR spectra suitable for DNN pure shift processing

Please use the dataset given in the **NMR/Kanamycin/3** along with the **Bruker** compatible pulse sequence **spin-echo\_pr** (a presaturation version); however, for the normal cases, where there is no significant water signal, **plw9** can be set to zero.